

Intravenous and Interstitial Photodynamic Laser Therapy Producing Groundbreaking Results in the Treatment of Cancer

Dr. med. Dipl. Chem. Michael Weber Weber Medical GmbH Sohnreystr. 4 37697 Lauenfoerde Germany info@webermedical.com / www.webermedical.com

+ 49 5273 367780





Content:

- 1. Basic Principles and Applications of Laser Therapy
- 2. Mechanisms of Photodynamic Therapy
- 3. Traditional Photosensitizers and their limitations
- 4. The new Nano ICG Photosensitizer
- 5. Applications and Protocol
- 6. First clinical results

What is Medical Laser Therapy?

- The therapy concept of LLLT has been developed in Germany and is based on the usage of **soft lasers**
- Unlike surgical lasers that destroy tissue and cells, soft lasers supply our body with energy to trigger numerous **regenerative processes**
- Medical laser therapy applies blue, green, yellow, ultraviolet, red and infrared lasers. Each color develops different effects which are evoked by the **stimulation of specific cellular components**, such as growth factors or cells of our immune system
- These effects have been proven by a large number of clinical studies







Main applications of Laser Therapy:

- Local Laser Therapy / Laserneedle Acupuncture for Pain Management and Tissue Regeneration
- Transcranial Laser Therapy
- Dermatological Therapy / Cosmetic Laser Therapy
- Interstitial / Intra-articular Laser Therapy
- Intravenous Laser Therapy
- Photodynamic Tumor Therapy

- In general, there are specific cellular structures that are able to absorb specific wavelengths (colors) of light (known as photoreceptors)
- The light stimulus gives a cellular signal affecting the chemical behavior, metabolism, movement and gene expression
- All associated enzymes and/or proteins are now affected
- This cascade event can ripple across an entire cell





Hamblin: http://photobiology.info/Hamblin.html



Absorption of Difference Light Wavelengths (Colors) in Mitochondria

- One example for the absorption of different colors within cells is the process in the mitochondrial respiratory chain [21]
- Complex 1 (NADH dehydrogenase) absorbs blue and ultraviolet light
- Complex 3 (cytochrome c reductase) absorbs green and yellow light
- Complex 4 (cytochrome c oxidase) absorbs red and infrared light





- Laser Light may help restore cellular homeostasis by inducing physiologically reparative activity for disease reversal in cancer and other degenerative diseases with minimal adverse side effects, and with potentially marked improvements in quality of life even in patients with advanced neoplasms
- Of major importance to achieve this is the activation and modulation of mitochondrial oxidative energy pathways
- Lasers have the potential to activate and modulate the production of ATP, GTP, AMPK and inositol pirophosphates P7-P8, not only through the respiratory chain but also through absorption and transportation of light by water
- A major goal for laser therapy in cancer is to control apoptosis (programmed cellular death) and differentiation, thus providing another cancer therapy tool

What is Photodynamic Therapy ?

- Photodynamic therapy is one of the most interesting and promising approaches in the treatment of various cancers
- The principle is the stimulation of a light sensitive drug which is injected into the blood (or given orally) and accumulates in cancer cells
- Tumor tissue is subsequently destroyed by irradiation with light of appropriate wavelength according to the absorption spectra of the various photosensitizers (by development of ROS)



What is Photodynamic Therapy ?

- Intravenous laser therapy in combination with photosensitive substances kills circulating tumor cells and stimulates the immune system
- Cancer tissue is targeted by local photodynamic therapy (i.e. interstitial or topical laser application)
- Common photosensitizers:
 - Liposomal Indocyanine Green (ICG)
 - Chlorin/Chlorophyllin
 - Curcumin (Turmeric)
 - Hypericin (St. John's Wort)
 - Riboflavin
 - Phycocyanin
 - Poly MVA
- Several chemo drugs are light-sensitive as well and can be enhanced by intravenous laser significantly (with lower dosage), i.e. Doxorubicin, Mitoxantron, Paclitaxel, Cisplatin or 5-FU



Light Absorption of Curcumin







Two individually non-toxic components brought together to cause harmful effects on cells and tissues:

- 1. Photosensitizing agent
- 2. Light of specific wavelength





















complement surviving components tumour cells PDT ○IL-1ß ○IL-8 neutrophils OIL-6 tumour cells PS necrosis apoptosis DC thrombus formation monocytes lymph node lymphocytes







Light Distribution and Cellular Response during PDT









Accumulation of PS in tumor

Result:

- Generation of reactive Oxygen Species (ROS)
- Apoptosis / Necrosis
- Elimination of PS from blood

Traditional (approved) Photosensitizers

- Hematoporphyrins, HpD
 - Derivatives of Hem (Photofrine and others)
- Chlorines
 - Derivatives of Chlorophyll
- Porphyrins
 - Synthetic Porphyrins

Platform	Drug	Substance	Manufacturer	Web site
Porphyrin	Photofrin®	HpD	Axcan Pharma, Inc.	www.axcan.com
Porphyrin	Levulan®	ALA	DUSA Pharmaceuticals, Inc.	www.dusapharma.com
Porphyrin	Metvix®	M-ALA	PhotoCure ASA	www.photocure.com
Porphyrin	Visudyne [®]	Vertiporfin	Novartis Pharmaceuticals	www.visudyne.com
Texaphyrin	Antrin®	Lutexaphyrin	Pharmacylics	www.pharmacyclics.com
Chlorin	Foscan®	Temoporfin	Biolitec Pharma Ltd.	www.bioletcpharma.com
Chlorin	LS11	Talaporfin	Light Science	www.lightsciences.com
Chlorin	Photochlor	HPPH	RPCI	www.roswellpark.org
Dye	Photosens [®]	Phthalocyanine	General Physics Institute	www.gpi.ru





Traditional (approved) Photosensitizers

TABLE 2. Clinically Applied Photosensitizers

PHOTOSENSITIZER	STRUCTURE	WAVELENGTH, nm	APPROVED	TRIALS	CANCER TYPES
Porfimer sodium (Photofrin) (HPD)	Porphyrin	630	Worldwide		Lung, esophagus, bile duct, bladder, brain, ovarian
ALA	Porphyrin precursor	635	Worldwide		Skin, bladder, brain, esophagus
ALA esters	Porphyrin precursor	635	Europe		Skin, bladder
Temoporfin (Foscan) (mTHPC)	Chlorine	652	Europe	United States	Head and neck, lung, brain, skin, bile duct
Verteporfin	Chlorine	690	Worldwide (AMD)	United Kingdom	Ophthalmic, pancreatic, skin
НРРН	Chlorin	665		United States	Head and neck, esophagus, lung
SnEt2 (Purlytin)	Chlorin	660		United States	Skin, breast
Talaporfin (LS11, MACE, NPe6)	Chlorin	660		United States	Liver, colon, brain
Ce6-PVP (Fotolon), Ce6 derivatives (Radachlorin, Photodithazine)	Chlorin	660		Belarus, Russia	Nasopharyngeal, sarcoma, brain
Silicon phthalocyanine (Pc4)	Phthalocyanine	675		United States	Cutaneous T-cell lymphoma
Padoporfin (TOOKAD)	Bacteriochlorin	762		United States	Prostate
Motexafin lutetium (Lutex)	Texaphyrin	732		United States	Breast



Traditional (approved) Photosensitizers

Photosensitizer	Type of diseases	Country			
(5-ALA)	Actinic keratosis,				
5-aminolevulinate	Basal cell carcinoma	U.S., EU			
Photofrin	Barrett's displasia	U.S., Canada, EU, UK			
Photofrin	Cervical cancer	Japan			
Photofrin	Endobronchial cancer	Canada, Most EU Countries, Japan, U.S.			
Photofrin	Esophageal cancer	Canada, Most EU Countries, Japan, U.S			
Photofrin	Gastric cancer	Japan			
Photofrin	Papillary bladder cancer	Canada			
Foscan	Head and neck cancer	EU, Norway, Iceland			
Verteporfin	Age-related Macular Degeneration	Canada, Most EU Countries, Japan, U.S.			



Limitations of traditional Photosensitizers







- Limited success by using red, blue and yellow lasers due to limited penetration depth (max 2 cm for red laser)
- Limited tumor size: max 2 cm
- Risk of burning and ulceration due to overdosage
- Light sensitivity
- Low success rate for liver and bone metastases
- No success in treatment of brain tumors



Key Features:

- First Photosensitizer in infrared range (infrared lasers have the best penetration depth into the body; even into bone tissue)
- Improved tumor targeting and accumulation with new nanotechnology
- Strong accumulation and tissue lifetime in tumors for several days: A series of laser treatments can be done after just one infusion/injection
- Large quantum yield of about 36 %
- Very safe substance with almost no side effects and negligible light sensitivity
- Highest quality standards guaranteed by a certified German sterile laboratory



Penetration of infrared laser through bones



What is Indocyanine Green (ICG)?

- Indocyanine Green (ICG) is a cyanine fluorescent dye used in medical diagnostics (for determining cardiac output, hepatic function, liver blood flow and for ophthalmic angiography)
- The hydrophilic, anionic tricarbocyanine photosensitizer has been developed by Kodak Research Labs and got FDA approval in 1956
- It has a peak spectral absorption at about 800 nm. These infrared frequencies penetrate retinal layers, allowing ICG angiography to image deeper patterns of circulation than fluorescein angiography
- ICG binds tightly to plasma proteins like albumin and is confined to the vascular system
- When ICG is administrated intravenously half-life is approx. 3 4 min
- Around 3,000 scientific papers on ICG have now been published worldwide



Fig 1: Absorption spectra of Indocyanine Green (ICG, 1,7-Bis[1,1-dimethyl-3-(4-sulfobutyl)-1H-benz[e]indol-2-yl]heptamethinium-betain sodium salt. A) Normalized spectra of a 10 mg/ml solution in water. B) Concentration dependence of the ICG spectra. Higher concentrations lead to a molecular assembly of the ICG molecules and a strong shift in the wavelength maximum could be observed. The effectivity of ICG decreased with the molecular assembly.





Limitations of regular ICG for Photodynamic Therapy:

- After intravenous injection, regular ICG binds to serum proteins almost completely (98 %) within 1-2 s; therefore, an accumulation in peripheral tissue, organs or even tumors is practically impossible (Cherrick et al., 1960)
- The tissue lifetime/stability of regular ICG in normal tissue and tumor tissue is about 10 min (Licha et al., 2002)
- The quantum yield of regular ICG in physiological NaCl and in plasma is about 3-5 % (Sutterer et al., 1996)
- → These properties do not allow any application in Photodynamic Therapy (PDT)



Solution: Nanoscale drug formulation

- One of the most useful advantages of nanoscale drug formulations is the utilization of the so-called enhanced permeability and retention effect (EPR)
- The concept of EPR describes all the mechanisms of passive transport of drug carriers of a defined size of about 200 nm or smaller, typically liposomes, micelles or nanoparticles to accumulate in tumor tissue much more than they do in healthy tissues



Why nanoscale drug formulations?

- Explanation: In order for tumor cells to grow quickly, they must stimulate the production of blood vessels by VEGF; additionally, the area surrounded by the fast-grown tumor is permanently inflamed; these newly formed tumor vessels are usually abnormal in form and architecture; they are poorly aligned defective endothelial cells with wide fenestrations; furthermore, tumor tissues usually lack effective lymphatic drainage



Development of Nano-ICG:

Liposomes and micelles were the first investigated nano systems for improved and better targeted introduction of photosensitizers into tumor cells. They have similar structures characterized by a core/shell model and thus can incorporate hydrophilic, hydrophobic and amphiphilic substances for delivery into cancer cells (Fig 3 and 4). The result is a considerable improved and selective uptake into tumor cells with an extended retention time of the photosensitizer.





The Transmission Electron Microscopy images visualizing the formulation confirm the PCS results of a monomodal and homogenous size distribution with an average hydrodynamic diameter below 100 nm.

Measurement with atomic force microscopy.













Outrinkes.

The new Nano ICG Photosensitizer



The baseline stability of Nano-ICG in normal muscle tissue is about 5 days. This is much longer than any other known formulation. The stability of Nano-ICG in tumorous tissue depends on the tumor type and tumor size; however latest measurements indicate it will be at least 3-4 days.





Accumulation in tumors: First measurements

Recent measurements now even confirm significant endogenous accumulation of Nano-ICG in liver carcinoma 3 days after infusion (interstitial measurement by infrared spectroscopy directly inside the liver). Thus, liver carcinoma and liver metastasis may effectively be treated (after infusion of 150 mg Nano-ICG) by local PDT using interstitial and external infrared irradiation.





- Besides local injection of Nano-ICG into the tumor it is now finally possible to achieve accumulation in the tumor by intravenous infusion of Nano-ICG
- After infusion and/or local injection the tumor is then irradiated by infrared laser light
- Irradiation of superficial tumors on or below the skin can be performed externally without penetration of the skin by using infrared laser sources with high power density
- Deeper tumors or metastases should be treated by interstitial laser technology
- Treatment of bladder, prostate or rectal cancers is done with special fiber optic catheters
- As infrared lasers are known to achieve good bone penetration this could also be a new treatment option for brain tumors and bone metastases
The Technology: Weberneedle® Endolaser with high power infrared lasers for interstitial PDT



MADE IN GERMANY



International Medical Approvals:

- CE Approval (Europe) •
- Health Canada •
- US-FDA for External Laser Therapy (PDT = off-label use) •
- TGA Australia for External Laser Therapy •
- Taiwan-FDA •

Development with Support of German Government and European Union















Application of Intravenous Photodynamic Therapy:



Demonstration video: <u>https://www.youtube.com/watch?v=Ika1GkefjSU&feature=youtu.be</u>



For application of intravenous lasers the physician can chose between a standardized 22G canula (blue) or our specific weberneedle® butterfly.

After placing the canula (or butterfly) into the vein the so-called laser catheter (sterile weberneedle® lasercath) is inserted into the canula.

The laser is then connected to the lasercath guiding the light directly into the flowing blood. No saline or other fluid is needed.





New Development: Y-Canula (3-way canula) for infusion and intravenous laser application







Application of interstitial Photodynamic Therapy:

Demonstration Video:

https://www.youtube.com/watch?v=dIBoH-EoRGE&list=PLwmaC9-LsMKDP89IjVZghqcmjpkIEfwau&index=11



After placing the needle in or close to the tumor the fiberoptic light guide is inserted into the needle.

The laser is then connected to the light guide. It directly targets the tumor without losing any energy in the skin barrier.

Applications of PDT with new Nano ICG





Interstitial PDT treatment with several infrared lasers (A, B)

PDT of breast cancer patient with interstitial and external laser

Applications of PDT with new Nano ICG









- Application of PDT in
- A) Tongue cancer
- B) Lymph metastases
- C) Nasopharyngeal Cancer
- D) Thyroid Cancer

Applications of PDT in Brain Tumors







Applications of PDT in Urology

Applications of PDT in Gastroenterology

Protocol for PDT with Nano ICG:

Dosage Information and Protocols:

Dosage/vial: 150mg/10ml

Application and 7-day treatment protocol:

Day 1: Dissolve 50 % of the vial (75 mg) in 500 ml glucose 5 % and infuse within 60 min.

Day 2: Repeat. If possible, do local injections into tumor*.

Day 2-7: Daily intravenous laser therapy** with red (50 mW), green (25 mW), blue (75 mW) and yellow (50 mW) laser as well as local irradiation of the tumor by interstitial and/or external infrared laser.

*If local injections of Nano-ICG into tumors are possible (i.e. breast cancer), injection of up to 30 mg (dissolved in 100 ml glucose 5 %) directly into the tumor followed by interstitial laser treatment is recommended. The remaining 120 mg (60 mg on day 1 and 60 mg on day 2) should still be given intravenously to generate high plasma concentration and long-lasting tumor accumulation.

** Add infrared to intravenous laser protocol on day 3.

6 cm right breast with lymph metastasis

to heal

Interstitial laser therapy after injection of diluted Nano ICG, 3 ml each needle)

Additional topical infrared irradiation and sonodynamic therapy

Result: Complete healing after 4 courses (every 3 months)

Ultrasound before therapy

Ultrasound after 4 treatment cycles

Case report: Breast cancer left side, diagnosis 11/2017 Photodynamic Therapy in May 2018

30.4.2018

6.5.2018

Patient: 71 years, USA

FINDINGS:

There are scattered fibroglandular elements in both breasts that could obscure a lesion on mammography. Images were obtained using Hologic Selenia Dimensions full field digital mammography with tomosynthesis (3D). Current study was also evaluated with a Computer Aided Detection (CAD) system. The patient reports a history of right breast outer hemisphere pain, both upper outer and lower outer quadrantst. Mammography demonstrates evidence of a spiculate mass/focal architectural distortion in the right upper outer quadrant. The margins are indistinct but the area appears to measure 3.5 cm in long axis. No associated suspicious appearing calcifications are identified.

There are separate bilateral benign appearing calcifications. No additional masses are identified.

First ultrasound 29.10.2018

Third ultrasound 06.11.2018 (after first treatment week)

2nd treatment course: 2.2.2019

Result of PET scan from 3.4.2019:

IMPRESSION:

1. Discrete focus of hypermetabolic activity in the region of the sigmoid colon without a discrete underlying CT abnormality. Although this finding may represent physiologic activity, this finding is concerning for an underlying colonic malignancy. Correlation with endoscopic evaluation is recommended.

2. With regards to patient's known right breast cancer, no evidence of hypermetabolic right breast mass or hypermetabolic metastatic disease.

Reviewed and Interpreted by: Navid Moradshahi 3/4/2019 2:40 PM

Breast cancer patient with ulceration and 4 treatment cycles (last one with Nano-ICG):

the to

First clinical results:

Testimonial (USA):

"We've had some exciting results so far:

One patient with multiple brain tumors along with metastases through liver, bones had a complete disappearance of all the tumors in her brain other than one that had shrunk substantially and all the cancer markers had greatly reduced indicating that the other tumors were dying off.

Also had a patient with tumors all over her lungs and that completely disappeared.

We've seen some great reductions in breast cancers as well where in one of the patients the large tumor completely disappeared."

First clinical results: Oesophageal Cancer

First clinical results: Oesophageal Cancer

to hear

First clinical results: Squamous Cell Carcinoma

Diagnosis: Base of Tongue squamous cell carcinoma, HPV positive with bi-lateral nodal involvement

"I was very unwell for the first 10 days on my return to the UK. I have no doubt this was a Herxheimer effect of so much dead cancer in my body. My urine was copper coloured for days but gradually cleared. I endured a terrible cough and couldn't sleep for several nights unless fully propped up in bed. I kept bringing up huge amounts of very brown coloured mucus that had a terrible smell. On 3 occasions large lumps of tumour came into my throat and I've attached photos with dates. The long tube (with my wife's finger in the photo for scale) was quite solid and about 6cm long!"

First clinical results: Squamous Cell Carcinoma

"It is now 3 weeks since we were with you and I no longer have any mucus, am hardly coughing and feel as if there's a lot of space in my throat. In some ways I feel that everything needs to settle down in me and I'm supporting my health with my diet embracing the Budwig protocol. I have read extensively on her scientific discoveries and approach to cancer. Nevertheless I am so pleased with what's happened to me with your care."

First clinical results: Lung Cancer

Case report from Malaysia:

"Lung Cancer patient did chemo 3 months before. 3 month later massive recurrent involving whole lung and pleura (stage 4). Surgeons gave no hone. So we tried new Nano

Surgeons gave no hope. So we tried new Nano ICG."

"I'm reporting from Malaysia [...] on the success case of my lung cancer patient (stage 4). [He is] one month after PDT with Nano ICG now near recovery."

Chest xray specialist said almost complete cure of lung and pleura.

BMG - Baja California Mexico (Jennifer Miele):

IMPORTANT INTEGRATION

Laser Therapy and PDT

- Integrating laser therapy and PDT into our clinical practice has shown impressive results. The synergy between therapies has noted increased responses and improvements.
- We continue our research and protocol integration for the results of the individual patient.

e 10

BMG - Baja California Mexico (Jennifer Miele):

"We have been working with the Nano ICG since late January 2019 with impressive and promising results ranging **from chronic infection**, lyme to late stage cancers."

64-year old female patient diagnosed with stage four Endometrial cancer. Metastatic disease with peritoneal tumor, liver tumor and right lung taken.

Report from Jan 2019: Liver tumor 5.2 CM by 5.6 CM Peritoneal tumor 2.2 CM X 2.8 CM March 2019: Liver tumor 2.0 CM X 2.4 CM Peritoneal tumor 2.7 CM X 0.6 CM

BMG - Baja California Mexico (Jennifer Miele):

19-year old Osteosarcoma patient with metastatic disease to liver.

PET CT scan reported SUV uptake of 14.2 in liver mass with 2.8 CM liver mass.

Follow up PET CT, after Nano ICG and PDT protocol, reported liver mass SUV 2.7 and mass size 1.7 X 1.2 CM.

BMG - Baja California Mexico (Jennifer Miele):

STAGE IV BREAST CANCER

57 y/o female, diagnosed in 2013 Diagnosis confirmed by biopsy – pathology report

PET CT Scan Diagnostic impression - September 2018 Metabolic tumor activity disseminated in liver and skeleton. Increased metabolism in Thyroid nodule. Metabolic activity in zone of the lung as well as mediastinal adenopathies.

Extreme pain due to metastatic activity in the lumbar/sacral area and in hips Could not walk without assistance of a walker or cane Abdominal swelling and discomfort Weakness, weight loss, lack of appetite, dizziness, migraines, abdominal distention

BMG - Baja California Mexico (Jennifer Miele):

- March 2018

CA 15-3 and CA 27-29 $\,$

reported at over 11,000

- July 2018

<u>CA15-3 =1586.8</u>

Ref. range less than 31.3 U/ml

<u>CA 27-29 =1589.8</u>

Ref range less than 38.6 U/ml

-September 2018

CA15-3 = 818.5

Ref. range less than 31.3 U/ml

<u>CA 27-29 = 742.5</u>

Ref range less than 38.6 U/ml

- November 2018 CA15-3 = 339.5Ref. range less than 31.3 U/ml CA 27-29 = 367.1Ref range less than 38.6 U/ml -December 2018 CA15-3 = 97.5Ref. range less than 31.3 U/ml CA 27-29 = 113.4Ref range less than 38.6 U/ml

BMG - Baja California Mexico (Jennifer Miele):

Comparative Abdominal ultrasound

December 2018

Metastatic liver nodules measuring 9-13-15-18 and 26mm maximum diameter.

January 2019

The classic metastatic nodules have disappeared. There are at present fibroid nodules in liver segment VII and VIII measuring 11-12 and 17 mm maximum diameter. The liver resembles a fibroid liver post-chemotherapy. There are no dilation of intra-nor extrahepatic biliary ducts observed

January 2019 - CT San of Chest

No noted metastatic or cancer activity is observed.


BMG - Baja California Mexico (Jennifer Miele):

After therapy:

- Liver enzymes normal range.
- CA markers decreased.
- C reactive protein returned to normal range.
- Noted improvement in image reports.
- Patient able to walk without assistance.
- No pain removal of prescribed pain medications.
- Weight stabilized
- Increased energy



BMG - Baja California Mexico (Jennifer Miele):

CHOLANGIOCARCINOMA

76 y/o female, diagnosed in 2017 Diagnosis confirmed by biopsy – pathology report

- Extreme abdominal pain and bloating
- Weakness
- Weight loss
- Lack of appetite
- Nausea
- Dizziness
- Jaundice
- Confusion
- Elevated bilirubin, liver enzymes, ammonia...

BMG - Baja California Mexico (Jennifer Miele):

CHOLANGIOCARCINOMA

November 2018

<u>CA19-9 = 11295.64</u>

Ref. range less than 37 U/ml

December 2018

<u>CA19-9 = 9003.21</u>

Ref. range less than 37 U/ml

January 2019

CA19-9 = 2742.43

Ref. range less than 37 U/ml

After therapy:

- Jaundice not present.
- Liver enzymes reduced.
- Bilirubin reducing to normal range.
- CA markers reducing.
- Appetite returned & weight stabilized.
- Patient notes increased energy
- No abdominal distention



First results in prostate cancer treatments:

- Until now very limited success with prostate cancer treatments
- 2019: Treatment with Nano-ICG and intraurethral and intrarectal irradiation
- \rightarrow 4 cases with promising results





BMG - Baja California Mexico (Jennifer Miele):

Prostate Cancer

65 Yr Old Male Patient

August 5 2019: PSA 17.31

August 7th 2019 Transrectal ultrasound 5.27 X 4.2 X 3.8 CM

Pereferica right 8 X 6 mm

Vascularized solid nodule

Sept 21 2019: PSA 7.03

Sept 9th 2019

4.2 X 4.3 X 3.64 CM

Pereferica right 5.8 x 4.6 mm

Reduced size of vascularized nodule



Patients Experience and Safety:



- PDT with Nano ICG is usually offered in a 1- or 2-week treatment course with daily treatment
- Safe application without significant side effects (however Photodynamic therapy can lead to local reactions such as pain and inflammation)
- One of the key advantages of Nano-ICG is that there is no increased light sensitivity (as opposed to photosensitizers like Chlorin or Hypericin)



- 1. Ailioaie et. al (2011): Laser Regeneration of Nerve Injuries in Children
- 2. Aluani, P. (2007): Intravenous Laser Therapy: Case Reports from the Field
- 3. Andrae, F. (2007): The Effects of intravascular Low Level Laser Therapy in the Scope of a Redifferentiation Therapy of Malignant Tumors
- 4. Bakeeva L, Manteiffel V, Rodichev E, Karu T. Formation of gigantic mitochondria in human blood lymphocytes under the effect of an He-Ne laser. Mol-Biol-Mosk. 1993 May-Jun; 27 (3):608-17
- 5. Boev S, Selivonenko V. The impact of the intravenous He-Ne-Laser therapy on the antioxidant system in patients with stable insertion angina and postinfarkt cardiosclerosis. Klin-Med-Mosk, 1997; 75,12:30-3
- 6. Brill G, Grigoriev S, Romanova T. Changes of leucocyte metabolism in He-Ne laser blood irradiation in vitro. Proceedings of SPIE. 1993; 1981:204-209
- 7. Dimitriev A, Iudin V, Aparov N, Matyrnov V. Effect of intravascular laser irradiation of the blood on blood cells in pancreatitis. Klin Med (Mosk.). 1989; 67, 5:108-110
- 8. Driianskaia V. The clinico-immunological effects of immunotherapy in patients with acute pyelonephritis. Lik-Sprava. 1997; Jul- Aug; (4):89-92
- 9. Dube A, Hausal H, Gupta PK. Modulation of macrophage structure and function by low level He-Ne irradiation. Photochem. Photobiol. Sci. 2003; 2, 8:851-855
- 10. Funk J, Kruse A, Kirchner H. Cytokine production after helium-neon laser irradiation in cultures of human peripheral blood mononuclear cells. Journal Photochem. Photobiol. Biology, 1992; 16, 3-4: 347-355
- 11. Gasparyan L. Laser Irradiation of the blood. Laser Partner Clinixperience All Volumes 2003:1-4
- 12. Gordon WC, Casey DM, Lukiw WJ, Bazan NG. DNA damage and repair in light-induced photoreceptor degeneration. Invest Ophthalmol Vis Sci. 2002;43(11):3511-21



- 13. Gulsoy M, Ozer G, Bozkulak O, Tabakoglu H, Aktas E, Deniz G, Ertan C. LLLT increases lymphocyte proliferation. Journal Photochem. Photobiol. Biology, 2006; 82, 3:199-202
- 14. Hamblin R, Viveiros J, Changming Y, Ahmadi A, Ganz R, Tolkoff J, helicobacter pylori accumulates photoactive porphyrins and is killed by visible light. Antimicrobial Agents and Chemotherapy, 2005; 49, 7:2822-2827
- 15. Heine H, Schaeg G: Origin and function of "roodlike structures" in mitochondria. Acta anat. 1979; 103:1-10
- 16. Heine H. Lehrbuch der biologischen Medizin. Stuttgart: Hippokrates, 3. Auflage 2007
- 17. Humpeler E, Mairbaurl H, Honigsmann. Effects of whole body UV-irradiation on oxygen delivery from the erythrocyte. Eur J Appl Physiol. 1982; 49:209-2014
- 18. Ivaniuta OM, Dziublik AIa, Skopichenko VN. Optimal effectiveness of complex treatment of patients with chronic obstructive bronchitis by intravascular laser irradiation of blood. Probl Tuberk. 1992;(5-6):21-4.
- 19. Karp G. Molekulare Zellbiologie. Heidelberg: Springer, 4. Auflage 2005
- 20. Karu T. Ten Lectures on Basic Science of Laser Phototherapy. Gangesber, Sweden: Prima Books AB (2007)
- 21. Karu T. The Science of Low-Power Laser Therapy. Amsterdam: Gordan and Breach Science Publishers, 1998
- 22. Karu, TI, Afanasyeva, NI. Cytochrome oxidase as primary photoacceptor for cultured cells in visible and new IR regions. Doklady Akad Nauk. 1995;342:693-695.
- 23. Kassak P, Sikurova L, Kvasnicka P, Bryszewska M. The response of Na/K-Atpase of human erythrocytes to green laser light treatment. Physiol Research, 2006; 55,2:189-194
- 24. Khotiaintsev K, Doger-Guerrero E, Glebova I, Svirid V, Sirenko J. Laser blood irradiation effect on electrophysiological characteristics of acute coronary syndrome patients. Proc. SPIE. 2929: 1996:132-137



- 25. Kipshidze N, Chapidze G, Bokhua M, Marsagishvili L. Effectiveness of blood irradiation using a Helium-Neon-Laser in the acute period of myocardial infarction. Sov-Med. 1990; 3:9-12
- 26. Kolarova H, Ditrichova D, Smolan S. Effect of He-Ne laser irradiation on phagocytotic activity of leukocytes in vitro. Acta-Univ-Palcki-Olmuc-Fac-Med. 1991; 129:127-132
- 27. Kozhura V, Dvoretskii S, Novoderzhkina I, Berezina T, Kirsanova A, Iakimento D, Kozinets GI. The effect of intravascular helium neon laser blood irradiation on the state of the compensatory processes in the acute period of hemorrhagic shock and after resuscitation. Anesteziol Reanimatol. 1993; 4:43-8
- 28. Kreisler M, Christoffers AB, Al-Haj H, Willershausen B, d'Hoedt B (2002). Low-level 809-nm diode laser-induced in vitro stimulation of the proliferation of human gingival fibroblasts. Lasers Surg Med 30(5):365–369
- 29. Ledin A, Dobkin V, Sadov A, Galichev K, Rzeutsky V. Soft-laser use in the preoperative preparation and postoperative treatment of patients with chronic lung abscesses. Proc. SPIE. 1999; 3829:2-5
- 30. Leonova G., Maistrovskaia O, Chudnovskii V. Helium-neon laser irradiation as inducer of interferon formation. Vopr-virosol. 1984; 39, 3:119-121
- 31. Lim WB, Kim JS, Ko YJ, Kwon H, Kim SW, Min HK, Kim O, Choi HR, Kim OJ. Effects of 635 nm light-emitting diode irradiation on angiogenesis in CoCl -exposed HUVECs. Lasers Surg Med. 2011;43(4):344-52.
- 32. Lindgard A, Hulten L, Svensson L, Soussi B. Irradiation at 634 nm releases nitric oxide from human monocytes. Lasers Med Sc 2007; 22:30-36
- 33. Litscher G./ Chiran/ Weber/ Ailioaie L./ Ailioaie C./Litscher D. (2013): Intravenous Laser Blood Irradiation Increases Efficacy of Etanercept in Selected Subtypes of Juvenile Idiopathic Arthritis - An Innovative Clinical Approach
- 34. Lutoshkin M, Tsypilev M, Lutoshkina M. Application of a Helium-Neon Laser (HNL) for the correction of renal function in patiets with chronic glomerulonephritis. Uro. Nefrol (Mosk.). 1993; 2:17-20
- 35. Manteifel V, Bakeeva L, Karu T. Ultrastructurel changes in chondriome of human lymphocytes after irradiation with He-Ne laser: appearance of giant mitochondria. Journal Photochem. Photobiol. Biology, 1997; 38, 25-30:



- 36. Manteifel V, Karu T. Structure of Mitochondria and Activity of their Rspiratory Chain in successive Generations of Yeast Cells exposed to He-Ne Laser Light. Biology Bulletin 2005; 32, 6:556-566
- 37. Meshalkin E. (ed.) Application of Direct Laser Irradiation in Experimental and Clinical Heart Surgery [in Russian], Novosibirsk: Nauka, 1981
- 38. Mi X, Chen J, Cen Y, Liang Z, Zhou L. A comparative study of 632,8 and 532 nm laser irradiation on some rheological factors in human blood in vitro. J. Photochem. Photobiol. B., 2004; 74,1:7-12
- 39. Miley G, Christensen. Ultraviolet blood irradiation therapy: Further studies in acute infections. American Journal of Surgery. 1947;73(4):486-493.
- 40. Miley, G. Effacacy of ultraviolet blood irradiation therapy and control of Staphylococcemias. American Journal of Surgery. 1942;64(3):313-322
- 41. Miley, G.: The Knott Technique of ultraviolet blood irradiation in acute pyogenic infections. New York State Journal of Medicine. 1942: 38- 46.Pathogen Deactivation
- 42. Mittermayr et al. (2007): Blue Laser Light Increases Perfusion of a Skin Flap Via Release of Nitric Oxide from Hemoglobin. Mol Med. 2007 Jan-Feb; 13 (1-2): 22-29:
- 43. Mouayed A, Fareed F, Ihsan F, Ahmad Y. Estimation of IgM & IgG values in the serum after intravenous irradiation of blood with diode laser. First UAE International Conference on Biological and Medical Physics, Al-Ain. 2005; Abstract No. 70
- 44. Noohi/ Javdani/ Kiavar (2008): Study of the Efficacy of Low Level Laser in Myocardial Perfusion in Patients with Chronic Stable Angina
- 45. Raggi, F./ Vallesi G. (2008): Intravenous Laser Blood Irradiation in Sports Medicine
- 46. Ramabhadran TV, F. T. (1976). In vivo induction of 4-thiouridine-cytidine adducts in tRNA of E. coli B/r by nearultraviolet radiation. Photochem Photobiol, 23(5), 315-21.
- 47. Ramabhadran TV, J. J. (1976). Mechanism of growth delay induced in Escherichia coli by near ultraviolet radiation. PNAS, 73(1), 59-63.



- 48. Safavi SM, Kazemi B, Esmaeili M, Fallah A, Modarresi A, Mir M. Effects of low-level He-Ne laser irradiation on the gene expression of IL-1beta, TNF-alpha, IFN-gamma, TGF-beta, bFGF, and PDGF in rat's gingiva. Lasers Med Sci. 2008 Jul;23(3):331-5. Epub 2007 Sep 5
- 49. Schmidt S, K. J. (2007). Process and Laboratory Scale UV Inactivation of Viruses and Bacteria Using an Innovative Coiled Tube Reactor. Chemical
- 50. Schulte-Uebbing, C. (2014): Intravenous Laser Therapy as an Integrative Approach in Immunology
- 51. Schulte-Uebbing, C. (2015): Intravenous Laser Therapy in an Integrative Therapy Concept for Treatment of Autoimmune Diseases (Hashimoto's Disease)
- 52. Schumm, N. (2006): Intravenous Laser Therapy: A Pilot Study for Treatment of Multiple Sclerosis
- 53. Schumm, N. (2008): Intravenous Laser Blood Irradiation in Multiple Sclerosis: A New Therapeutic Procedure with Significant Improvement in Quality of Life
- 54. Siniukhin VN, Ianenko EK, Safanov RM, Khamaganova EG, Borisik VI. The effect of endovascular helium-neon laser therapy on the immune status of patients with acute calculous pyelonephritis. Urol Nefrol (Mosk). 1996 Nov-Dec;(6):9-11.
- 55. Siposan G, Lukacs A. Relative variation to received dose of some erythrocytic and leukocytic indices of huma blood as a result of low-level-laser radiation. Journal of Clinical Laser Medicine & Surgery. 2001; 19, 2:89-103
- 56. SpasovA, Nedogoda V, Konan K, Kucheriavenko A. Effect of the intravenous laser blood irradiation on efficacy of drug preparations. Eksp Klin Farmakol. 2000; 63, 5:65-7
- 57. Stadler I, Evans R, Kolb B, Naim J, Narayan V, Buehner N, Lazafame R. In vitro effects of low level laser irradiation at 660 nm on peripheral blood lymphocytes. Lasers Surg. Med. 2000; 27, 3:255-261
- 58. Stroev E, LarionovV, Grigoreva L, Makarova V, Dubinina I. The treatment of diabetic angiopathies by endovascular low-intensity laser irradiation. Probl-Endokrinol-Mosk. 1990; 36, 6:23-5



- 59. Tuner J, Hode L. Laser Therapy Clinical Practice and Scientific Background. Grängesberg: Prima Books AB, 2002
- 60. Vasa/Breitschopf/Zeiher/Dimmeler: Nitric Oxide Activates Telomerase and Delays Endothelial Cell Senescence
- 61. Vinck E, Cagnie B, Cornelissen M, Declerque H, Cambier D. Green light emitting diode Irridation enhances Fibroblast Growth impaired by high glucose levels. Photomedicine and laser surgery. 2005, 23, 2:167-171
- 62. Weber et al. (2010): Intravenous Laser Blood Irradiation: Introduction of a New Therapy
- 63. Wieden, T. (2009): Fibromyalgia in Pain Therapy. Mechanisms and Treatment Opportunities in Laser Therapy
- 64. Wirz, A. (2008): Frequencies With Consequences. Pilot Study for Treating Diabetes
- 65. Wirz/ Baumgartner/ Burger/ Gerber (2008): Intravenous Laser Therapy in the Treatment of Horses. Results of a Multi-Center Pilot Study
- 66. Xuejuan G, Xing D. Molecular mechanisms of cell proliferation induced by low power laser irradiation. J Biomed Sci. 2009;16(1):4
- 67. Yang Z , Wu Y, Zhang H, Jin P, Wang W, Hou J, Wei Y, Hu S. Low-level laser irradiation alters cardiac cytokine expression following acute myocardial infarction: a potential mechanism for laser therapy. Photomed Laser Surg. 2011 Jun;29(6):391-8. doi: 10.1089/pho.2010.2866. Epub 2011 Feb 24
- Yu H, Chang K, Yu C, Chen J, Chen G. Low-energy helium-neon laser irradiation stimulates interleukin-1 alpha and interleukin-8 release from cultured human keratinocytes. Journal of investigative Dermatology. 1996; 107:593-596
- 69. Zhang H, Hou JF, Shen Y, Wang W, Wei YJ, Hu S., Low level laser irradiation precondition to create friendly milieu of infarcted myocardium and enhance early survival of transplanted bone marrow cells. J Cell Mol Med. 2010 Jul;14(7):1975-87. doi: 10.1111/j.1582-4934.2009.00886.x. Epub 2009 Sep 1
- 70. Zuern, I. (2016): Pilot Study on Treatment of Chronic Lyme Disease with Yellow and Blue Laser

