
Handy Cure s'®

Treatment Area

- Where possible, treatment should be provided over bare skin.
- During treatment, the device should be held in a hovering state over the area where the pain is felt.
- If the treatment area is wider, the device should be moved (scanned) around the treatment area
- If required, treatment can be provided over clothing or over wound dressing.

experience relief.



	Treatment area and number of treatments	Program week 1	Program weeks 2 and 3
1	Achilles' tendon - 2 scanning sessions	Variable	50 Hz
2	Ankle - 2 sessions, front and back	Variable	50 Hz
3	Elbow joint (Tennis Elbow; Golf Elbow) - 2 sessions, front and back	Variable	50 Hz
4	Fingers (Non-deforming arthritis; post-traumatic processes) - 2 sessions, front and back	Variable	50 Hz
5	Foot, toes - 3 sessions, both sides of ankle and scanning top of foot	Variable	50 Hz
6	Heel (Calcaneal spur) - 2 sessions, back of ankle and under heel	Variable	50 Hz
7	Lower Back Pain (disci intervertebrales problems, Ischialgia) - 2 sessions scanning lower back, 1 session scanning back of thigh	Variable	50 Hz
8	Migraines, headaches - 2 sessions scanning back of neck	Variable	50 Hz
9	Shin - 2 scanning sessions	Variable	50 Hz
10	Spine (Vertebral Column Osteochondrosis) - 4 sessions scanning affected area	Variable	50 Hz
11	Wrist (Writing Syndrome, Radiocarpal Articulation)-2 sessions, front and back	Variable	50 Hz
12	Knee - 4 sessions around knee	Variable	50 Hz
13	Hip joint - 4 scanning sessions	Variable	5 Hz
14	Shoulder - 2 sessions scanning front and 2 sessions scanning back of shoulder	Variable	5 Hz
*	Arthritis, ligament or tendon injury, cartilage abrasion – treatment in the pain area as indicated in the table. For chronic pains, including; muscle, tendon pain and other injuries not mentioned in the above table, we recommend performing 2 treatment sessions with the Variable program for two weeks or until you		

Handy Cure s'®



