

Local PBM Applications

From now on, a **separate controller** can operate up to two applicators simultaneously


Mouth/Throat Applicator UV/Blue

Diode type: LED Wavelengths:

- 14x 447 nm (Blue)
- 14x 370 nm (UV)

Max. output: 3.4 W





- 3x 650 nm (Red)
- 3x 532 nm (Green)
- 3x 450 nm (Blue)
- 3x 589 nm (Yellow)

Max. output: 1.4 W



Nasal probe - Red

Diode type: Laser Wavelength: 2x 650 nm Max. output: 135 mW



Diode type: Laser

Wavelength: 2x 650 nm Max. output: 120 mW



Diode type: LED Wavelengths:

- 1x 447 nm (Blue)
- 1x 370 nm (UV)

Max. output: 210 mW



Controller

Select duration and intensity at the press of a button.

Duration: 5 - 10 - 15 - 20 - 25 - 30 - 35 - 40 - 45 - 50 - 55 - 60 minutes

Intensity: 3 levels

For more info and orders: info@weberlaser.com www.weberlaser.com

Quick Guide

Controller

1. Charge the controller with the included USB cable and power plug.



The charging port is located on the bottom of the device (1).

The charging process is completed when the bars of the battery symbol stop blinking.

Your controller powers the different applicators.

It has two ports so that you can use two applicators simultaneously.



- 2. Plug the desired applicator(s) into the ports provided for this purpose (2) and apply them to the desired body areas (e.g. into the nose/ears/mouth, onto the skin, etc.).
- 3. Press the **On/Off button** for 3 seconds to turn the controller on.
- 4. Press the **time button (T)** to set the desired duration: 5 10 15 20 25 30 35 40 45 50 55 60 minutes.
- 5. Press the **power button (P)** to start the session.
 - During the session you have the option to press the power button (P) any time in order to adjust the intensity (3 levels).
- 6. You can end the session early by pressing the On/Off button (hold for 3 seconds).



Recommendations for use

We recommend 1-2 sessions/day with 20-30 minutes per session.

At the beginning, start with low intensity (Level 1) and/or shorter duration (5-10 minutes) and observe your reaction. You can then gradually increase intensity and duration from session to session depending on how you feel.

If you require any help with the set-up, use or maintenance or to report unexpected operation or events, please contact the general service address: info@wmedicalsystems.com