

User Manual

To Enjoy A Healthy Life



Instruction Manual

Congratulations on taking control of your pain !
Now you can use the LED Red Light Therapy to provide direct, non-addictive, drug free pain relief.

Efficient Cells Can Enhance Health

Please read the instructional manual to get the most out of your device.

Features:

1. Perfect performance for body treatments like weight loss and pain relief.
2. Portable design so you can use it in different situations.
3. LEDs with wavelength 660nm 850nm infrared.
4. We provide 1 year warranty for belt.




Weight loss



Pain relief

Parts Of The Light Therapy Belt



105pcs
Light Therapy Belt



300pcs
Light Therapy Belt

Plugs into DC Power Port on Belt

Remote Control

Display Screen



Remote Control and Cable

Plugs into DC Port on Remote Control Cable



Plugs into USB Power Source

*DC USB Power Adapter Cable

Only 105pcs
Red Light Therapy Belt

Wall Plug (AU,EU,UK,US)



Plugs into DC Port on Remote Control Cable

AC Power Adapter

Operating Instructions

Please read the operating instructions before using the Light Therapy Belt. Keep this manual for reference.

•Powering The Light Therapy Belt With The Wall Adapter

1. Connect the male end of the remote control cable to the DC power port exiting the light therapy belt.
2. Connect the wall power adapter to the female end of the remote control cable.
3. Plug the power adapter into a standard wall outlet.

•Powering The Light Therapy Belt With The USB Adapter

1. Connect the male end of the remote control cable to the DC power port exiting the light therapy belt.
2. Connect the USB adapter to the female end of the remote control cable.
3. Insert the other end of the USB adapter into a USB device, such as a laptop or computer. Alternately, you may insert the USB connector into a standard mobile phone 5V-3A or 5V-4A USB charging adapter and then plug the adapter into a wall outlet.

WARNING: The light therapy belt cannot be used with a 12V USB charging adapter.

•Powering The Light Therapy Belt With An External Battery

The light therapy belt supports the use of a 5V external battery, also called a power bank. The external battery is not included and must be purchased separately from other retailers. To power the light therapy belt with an external battery:

1. Connect the male end of the remote control cable to the DC power port exiting the light therapy belt.
2. Connect the external battery pack to the female end of the remote control cable, via the cable included with the battery pack.
3. Place the battery in the pocket at the front of the belt.

•Using The Control Panel

Mode Setting:
660nm/850nm/ALL

Brightness Setting:
L0,L1,L2,L3,L4

Timer Setting:
5,10,15,20,25,30 Mins

Switch:
Open/10Hz Flicker/OFF



105pcs

Red Light Therapy Belt

Increase Brightness
Setting:L0,L1,L2,L3,L4

Timer Setting:
5,10,15,20,25,30 Mins

Decrease Brightness
Setting:L0,L1,L2,L3,L4

Switch:
Open/10Hz Flicker/OFF



300pcs

Red Light Therapy Belt

•Intensity Options

Five different intensity settings allow you to choose the right setting for your particular need. Each of the five settings offers a different power (in Watts) to help reduce inflammation, increase blood flow, relieve pain, and lose fat.

| INTENSITY LEVEL | 105Pcs POWER (WATTS) | 300Pcs POWER (WATTS) |
|-----------------|----------------------|----------------------|
| L0 | 0.9W | 9W |
| L1 | 4.9W | 17W |
| L2 | 8.5W | 25.7W |
| L3 | 12.3W | 34.4W |
| L4 | 16.0W | 40.9W |

Support two ways of usage

The red light belt can be directly plugged into the power supply or plugged into the power bank for use.

1 DC Charging

Input Voltage: AC 100-240V, 50/60Hz



105pcs or 300pcs
Red Light Therapy Belt

2 Power Bank Cord (*Only 105pcs Red Light Therapy Belt are Supported)

Input Voltage: DC 12V. With Power Bank Cord, you can connect with any power bank (not included) which is DC 12V, and put the power bank to the bag.



Only 105pcs
Red Light Therapy Belt

Product Specification

| | | |
|------------------------------------|--|------------------------|
| LED Number: | 105Pcs | 300Pcs |
| Product Size: | 49.6*6.8 inch | 63.7*11.8 inch |
| Net Weight: | 2.2 lbs | 3.3 lbs |
| Frame Color: | Black/Blue | Black |
| Wavelength: | RED/NIR | RED/NIR |
| Function: | Body slim/ Pain Relief | Body slim/ Pain Relief |
| Customize Service: | Logo/Shape/Wavelength | Logo/Shape/Wavelength |
| AC Power Rated Input: | 100-240 V ~50/60 Hz | 100-240 V ~50/60 Hz |
| DC Power Rated Output: | 12V 2A | 12V 4A |
| DC USB Power Adapter Cable Input: | 5V 3A or 5V 4A | Not support |
| DC USB Power Adapter Cable Output: | 12V 1.67A | Not support |
| *Power Bank Output: | 5V 3A or 5V 4A (*Not included; sold separately from other retailers.) | Not support |

Caution

*Don't use the light around water

*Don't stare directly at the light use protective glasses while using the device

*Don't treat near the eyes

Equipment Maintenance

Cleaning And Disinfecting:

Use only a soft dry cloth to clean the exterior. Never use an abrasive cleanser or submerge in water or any other liquid. This device cannot be sterilized by any liquid or autoclave method.

Storage: Store the unit in a cool, dry place, away from direct light and extreme heat.

Tips For Using The Light Therapy Belt

•Using The Light Therapy Belt

- Place the light therapy belt in position on your waist before powering on the device.
- Turn off the light therapy belt before removing it from your waist.
- Do not use continuously for more than 30 minutes, and wait at least 2 hours between 30-minute sessions.
- Use the light therapy belt daily for optimal results.
- We recommend washing and drying your hands thoroughly and only using over clean clothing and skin to keep the therapy surface clean.
- While the light therapy belt is primarily designed for use on the waist and lower back, it may also be used on the upper back, shoulders, arms, legs, hips, or feet. Do not use around the head or neck.
- Do not wear accessories, such as belts or jewelry, under the light therapy belt.
- The remote control cable and power adapter or USB adapter should be disconnected from the belt for storage or travel.

•Important Safety Information

- Never use the light therapy belt under a blanket or pillow. Excessive heat can potentially cause fire, electrocution, or injury to people or property.
- Keep the light therapy belt , remote control, and power cables away from heat sources and fire.
- The device is not waterproof. Keep away from liquids. Do not immerse in water.
- Never use the device or plug/unplug the power cables with wet hands.
- Always power off the device after each use.
- Do not use or wrap the light therapy belt or cables around your head or neck.

Note:

Indoor use only: Light is not waterproof, can not used in water or in wet places.

Repairing: If your light need after-sale service, please contact us directly

Installation : To avoid the danger of electric shock, please turn off the power switch before installation.

